

# DRE Milano T50

## Power Procedure Chair

### Fully electric, multi-function power procedure chair

The DRE Milano T50 is a fully electric, multi-function chair that compares with popular brands in terms of style and functionality, at a fraction of the cost. With four actuators to adjust height, backrest, leg rest and Trendelenburg positioning, clinicians are ensured a clear path to the patient. The chair's low height capabilities give the patient easy access in and out of the chair.

The padded, anatomically-configured mattress foam and adjustable or changeable headrest make the Milano T50 Power Procedure Chair naturally adaptable to a variety of lengthy office-based procedures. It also comes standard with hand and foot controls with smart return, allowing clinicians to easily maneuver the chair into a neutral starting position.

#### FEATURES

- Fully electric operation with four actuators to adjust height, backrest, leg rest and Trendelenburg positioning.
- Upholstery is designed for medical use and is resistant to water, air and disinfectants.
- Standard handheld remote control and foot pedal include automatic reset function.
- Anatomic mattress is made of top-quality foam and includes a five-year warranty.
- Adjustable and removable headrest.
- Constructed using latex-free, non-toxic and non-carcinogenic materials.
- Seat and back rails fit a variety of medical accessories.
- Magnetic back remote control easily attaches to the side of the chair.
- ABS coated steel frame gives the chair a smooth finish while protecting the inner components.



*Shown with optional medical arm-rests*

390-2019-10-28

#### SPECIFICATIONS



Chair Weight: 196 lbs. (89kg)

Maximum Capacity Load:  
440 lbs. (200kg)



Adjustable Height:  
24" (62cm) to 36" (92cm)

Chair Dimensions:  
74" x 24" x 34"  
(190 x 62 x 86cm)



Trendelenburg: 25°

Backrest: 75°

Leg rest: 65°



Voltage:  
110 and 220 V cc

#### COLOR OPTIONS



Blue



Cream



Chrome



Black

**WARRANTY:** 2 Years

## MULTI-SPECIALTY USE

- › General Medical
- › Dermatology
- › Podiatry
- › Plastic, Cosmetic and Aesthetics, OB/GYN
- › ENT/Otolaryngology
- › Oral & Maxillofacial
- › Dental
- › Orthodontics

*Shown with optional accessories*



## OPTIONAL ACCESSORIES

- › I.V. Rod
- › Curtain-holder
- › Medical Arm-est
- › Articulating Arm Support
- › Jointed Thigh Support
- › Hand-set Control
- › Fixed Clamp
- › Rotating Clamp
- › Colposcope Support
- › Narrow Headrest
- › Articulating Knee Crutch
- › Castors
- › Instrument Tray

