

DRE Premio M3

Manual Hospital Bed

The DRE Premio M3 is a manual hospital bed with three movements: back rest, leg rest, and hi-low. Each table movement is easily controlled by cranks. The M3 features removable head, foot, and bed boards that are easy to clean. Includes collapsible railing and four universal casters with individual breaks.

FEATURES

- › Steel framework, top class ABS & PP plastic
- › Three functions controlled by cranks: Backrest, Leg Rest, and Hi-Low
- › Head, foot, and bed board are detachable and easy to clean
- › 6 column, collapsible aluminum alloy protective railing
- › 4 universal casters with individual breaks
- › Standard accessory: IV Pole



SPECIFICATIONS



Net Weight:
225 lbs (102kg)



Overall Size:
85" x 37.8" x 22.6 - 29.5"

Bed Board Dimensions:
75.2" x 32.7"

Hi-Low Adjustment:
29.5" - 22.6"



Backrest:
0-65° ± 5°

Leg Rest:
0-45° ± 5°



2358-2020-06-02