

Avante SLS 2500

LED Surgical Light

Provides a shadowless, natural color temperature

The Avante SLS 2500 LED Mobile Lighting System is a compact system. The shadowless lamps provide a color temperature of 3,500 k ~ 5,000 k, reflecting the actual colors of human tissue.

FEATURES

- Dimming knob located on the side of the lighthouse for quick and easy adjustments.
- Multipath light ensures that a single damaged LED in one path will not affect operation.
- The SLS 2500 is available as a ceiling mount! Ask your representative for details.
- Universal appending system provides maximal adjusting range, suitable for any operating room. Appending arm is made of nova alloy, with light-strong structure, convenient operation and accurate positioning.
- Lamp cap is totally enclosed with a streamlined design, complying with requirements of sterilizing and laminar flow cleaning in all kinds of operating rooms.



SPECIFICATIONS



Illumination Intensity:
60,000 lux

Color Temperature:
3,500 ~ 5,000 Kelvin

CRI - Color Rendering Index:
85 - 95

Diameter of Light Spot:
7"

Illumination Depth:
23.5"

Power:
20 W



Light Options:



Wall-Mounted



Mobile



Ceiling Mount

471-2019-08-26