

DRE SLS 9000

LED Surgical Light

Efficient surgery lighting system provides extremely cool output

The SLS 9000 is a shadowless light with a 50,000 hour rated LED life. Available in multiple configurations to fit your operating room needs.

FEATURES

- Cost efficient surgery lighting system.
- Available in a single or dual head ceiling mount, single wall mount or mobile floor model.
- Cool output LED light source.
- LCD touch screen operation measures illumination intensity, color temperature and Chromogenic index.



LCD touch screen operation



SPECIFICATIONS



Spot Diameter:
120 - 350mm

Color Temperature:
3500°-5000° Kelvin

LED Life:
50,000-hour

Illuminance:
90,000 lux
(9,000 foot candles)
per light head



Available in ceiling, floor, and wall configurations

482-2019-07-30