

DRE SLS 9000

Portable Floor Light

Efficient surgery lighting system provides extremely cool output

The SLS 9000 Portable LED light is energy saving and eco-friendly, with an ultra-long service life, up to 30,000 to 50,000 hours.

The CAD designed multi-LED-lighting-beam produces lighting depth more than 55 inches, with up to 90,000 lux of central illumination. New alloy construction and a simple suspension system is easy to control and be positioned accurately. This easily adjustable light system is suitable for pendant, wall-mounted and floor-standing lamps, meeting requirements of multiple operating rooms.



SPECIFICATIONS



Diameter:
20"

Color Temperature:
3500°-5000° Kelvin

LED Life:
50,000-hour

Illuminance:
90,000 lux
(9,000 foot candles)
per light head

Color Rendering Index:
93%

Diameter of Light Spot:
5" - 14"

LED Power:
55W

Illumination Depth:
28" - 59"

LED Service Life:
50,000 hours



*Easy-to-use
control panel*

897-2019-08-14