

DRE Waveline EZ MAX


Vital Signs Monitor

High resolution color touchscreen

The DRE Waveline EZ MAX is a full-featured touchscreen monitor designed for a variety of settings. The customizable menu allows operators to easily navigate between the parameters of the Waveline EZ MAX, including 3/5-lead ECG, RESP, NIBP, SpO₂, 2-TEMP and PR. The DRE Waveline EZ MAX features a no-fan design, ideal for quiet recovery and sterile environments.

FEATURES

- › 3/5-lead ECG, RESP, NIBP, SpO₂, 2-TEMP, PR
- › 10.4" touchscreen color TFT-LCD monitor
- › No-fan design
- › Arrhythmia, S-T Analysis, OxyCRG, Large Font Display
- › USB data storage and review
- › Bi-directional Communications with Central Monitoring System
- › Waveform color adjustable
- › Customizable alarm setting
- › Optional Parameters: Nellcor OxiMax™ SpO₂, 2-IBP, ICP, C.O. Respiration CO₂ (Mainstream and Sidestream)



SPECIFICATIONS



Height:
10" (254mm)

Width:
6.3" (161mm)

Length:
12" (303mm)

Display:
10.4" touchscreen color TFT,
800 x 600 resolution



Power:
100-240V,
50/60Hz

914-2019-10-25