

Hill Rom Beds for Safety, Quality and Comfort. Well-Designed, Clinically Effective and Easy to Use.

Specifications

- 4-motor bed
- Hi-Lo travel time: 21 seconds
- Sleep surface: 36" x 82"
- Overall length (without bumpers): 36" x 91"
- Bed frame clearance with 5" casters: 7.5"
- Low position: 17"
- High position: 30"
- Central brake and steer
- Automatic night light
- Instant CPR release



Features

- All beds meet or exceed manufacturer's specifications.
- We provide JCAHO required documentation.
- Paint color of choice (light or dark).
- Fully function tested.
- Basic bed function: Head, Knee, Hi-low.
- Orthopedic / Hardpan surface.
- Side rails include new panels and switches.
- New head and foot boards, your choice of laminate colors.
- New motor couplings and capacitors.
- Personalization with hospital logo upon request.
- 1 year warranty.



Hill Rom Advance Series Hybrid Hospital Bed

Equipment for the way *you* operate

Mattress Specifications

- All new and customized to meet the following:
 - » Dimensions: 84" L x 36" W x 5" D (for Hill Rom 1105)
 - » Dark blue water resistant nylon top half
 - » Heavy-duty gray vinyl bottom half
 - » Fire retardant to Cal 129 specifications



All new mattress!