

Refurbished Hill-Rom TotalCare hospital beds have advanced features that provide maximum comfort to patients and help reduce strain on nurses and physicians.

Perfect for hospitals, surgery centers and teaching facilities!

Product Details

- Professionally refurbished to the original manufacturer's specifications
- Total length: 92.5"
- 5" casters
- Bed lift capacity: 500 lbs.
- Trendelenburg/Reverse Trendelenburg
- Central brake and steer
- Sleep surface: 35" (w) x 83" (l)
- Height range: 14.5" to 33.5"
- Full chair position
- JCAHO-required documentation

Options

- Side guard communications
- Percussion and vibration module
- Rotation module
- Two-layer therapeutic surface
- Service contracts for preventative maintenance

